

MINI MASSAGER

Instruction Manual

8081440

INCLUDED

- Mini massager
- Instruction manual

BATTERY INSTALLATION

1. Turn the top of the massager counter-clockwise to remove the cover.
2. Install 3 AAA batteries (sold separately), aligned to the correct polarities (+/-).
3. Replace the cover according to the markings on the side of the product for placement.

OPERATING INSTRUCTIONS

1. Place the massager on the part of the body that you'd like to massage.
2. Press the centre button to turn vibration ON. Press again to turn OFF.

TROUBLESHOOTING

The device won't turn ON or it isn't vibrating powerfully:

- Replace the batteries
- Ensure the batteries are inserted properly

SPECIFICATIONS

Battery:	3 AAA (sold separately)
Size:	89 x 89 x 90 mm
Weight:	148g
Material:	ABS & silicone
Power Output:	1.3VA

CARE AND MAINTENANCE

- Use a soft cloth dampened with alcohol to clean the body
- Keep away from water, rodents and extreme temperatures
- **DO NOT** place heavy objects on top of this device
- **DO NOT** attempt to repair or alter the device on your own, please ask a professional to repair the device



THE SOURCE - WARRANTY

The Source warrants that this product will be free from defects in materials and workmanship for a period of 6 months from the date of purchase. Within this period, simply take the product and your proof of purchase to any The Source store or participating dealer and the product will be replaced (where available) without charge. Any product which has been subject to misuse or accidental damage is excluded from this warranty.

This warranty is only applicable to a product purchased through The Source company-owned stores or participating dealers in Canada where the warranty is included with the product. While this warranty does not confer any legal rights other than those set out above, you may have additional statutory rights which will vary under the laws of the various countries, states, provinces and other governmental entities in which The Source operates. This warranty is subject to all statutory rights you may have in Canada.

Imported by:
The Source, Barrie, Ontario, Canada, L4M 4W5
Manufactured in China
www.thesource.ca
© 2019 The Source. All rights reserved.

